

Warm.A.Pet™

In Car Heat Pad Controller

To adjust the heat of a Warm.A.Pet Small Deluxe, Large Deluxe or Joey Pouch heat pad when used in conjunction with the matching Warm.A.Pet In-Car Lead

Instructions for Use

- 1/. Attach the appropriate in-car lead to the heat pad.
(Replacing the normal source of power)
- 2/. Plug the cigarette lighter end of the in-car lead into the socket of the controller.
- 3/. Plug the cigarette lighter end of the Controller unit into the cigarette lighter facility in your car.
- 4/. Turn the control knob of the Controller unit, initially, to the high position (**H** position).
- 5/. Place the heat pad into the animals' bedding and (always) place a thermometer probe underneath where the animal will be, in the bedding to monitor the level of warmth generated by the heat pad on this initial (high) setting.
- 6/. Closely monitor the temperature and if/when required turn the control unit knob anti-clockwise towards low (**L**) to reduce the level of warmth from the heat pad.

Note: If the control unit is left set on high (**H**) then the maximum temperature will be attained which may be too hot for the animal.

Manufactured in Australia

by

Homark Pty Ltd

ABN: 23 057 128 300

62 Berkeley Street

Huntingdale Vic 3166

03 9544 1471 Fax 03 9544 6799